

DATA SHEET

Revision No.: ZJ0006 Issue Date: Oct 1st, 2011

Campylobacter Jejuni Real Time PCR Kit User Manual For In Vitro Diagnostic Use Only

REF DD-0040-01

For use with LightCycler1.0/2.0 Instrument



Gentaur Molecular Products Voortstraat 40 1910 Kampenhout, Belgium

1. Intended Use

C.jejuni real time PCR Kit is used for the detection of C.jejuni in stool or water samples by LightCycler (Roche) real time PCR systems.

2. Principle of Real-Time PCR

The principle of the real-time detection is based on the fluorogenic 5'nuclease assay. During the PCR reaction, the DNA polymerase cleaves the probe at the 5' end and separates the reporter dye from the quencher dye only when the probe hybridizes to the target DNA. This cleavage results in the fluorescent signal generated by the cleaved reporter dye, which is monitored real-time by the PCR detection system. The PCR cycle at which an increase in the fluorescence signal is detected initially (Ct) is proportional to the amount of the specific PCR product. Monitoring the fluorescence intensities during Real Time allows the detection of the accumulating product without having to re-open the reaction tube after the amplification.

3. Product Description

Campylobacter jejuni is a Gram-negative slender, curved, and motile rod. It is the leading cause of food poisoning, being three times more common than Salmonella. This bacterium is now recognized as an important enteric pathogen. Usually C.jejuni outbreaks are small (less than 50 people). C.jejuni infection causes diarrhea, which may be watery or sticky and can contain blood and fecal leukocytes. Other symptoms often present are fever, abdominal pain, nausea, headache and muscle pain

C.jejuni real time PCR kit contains a specific ready-to-use system for the detection of C.jejuni by polymerase chain reaction (PCR) in the real-time PCR system. The master contains reagents and enzymes for the specific amplification of C.jejuni DNA. Fluorescence is emitted and measured by the real time systems' optical unit during PCR. The detection of amplified C.jejuni DNA fragment is performed in fluorimeter channel 530nm with the fluorescent quencher BHQ1. DNA extraction buffer is available in the kit. In addition, the kit contains a system to identify possible PCR inhibition by measuring the 560nm fluorescence of the internal control (IC). An external positive control (1×10⁷ copies/ml) contained, allows the determination of the gene load.

4. Kit Contents

Ref.	Type of reagent	Presentation 25rxns
1	DNA extraction buffer	2 vials, 1.5ml
2	C.jejuni Reaction Mix	1 vial, 450μl
3	PCR Enzyme Mix	1 vial, 12μl
4	Molecular Grade Water	1 vial, 400μl
5	Internal Control (IC)	1 vial, 30µl
6	C.jejuni Positive control(1×10 ⁷ Copies/ml)	1 vial, 30µl

Analysis sensitivity: 1×10⁴ copies/ml

5. Storage

- All reagents should be stored at -20°C. Storage at +4°C is not recommended.
- All reagents can be used until the expiration date indicated on the kit label.
 Repeated thawing and freezing (>3x) should be avoided, as this may reduce the sensitivity of
- Cool all reagents during the working steps
- Reaction mix should be stored in the dark

6. Additionally Required Materials and Devices

- Biological cabinet
 Real time PCR system
- Desktop microcentrifuge for "eppendorf" type tubes (RCF max. 16,000 x g)
- Vortex mixer
- Real time PCR reaction tubes/plates
- Cryo-container
- Pipets $(0.5\mu l 1000\mu l)$
- · Sterile filter tips for micro pipets
- Sterile microtubes
- Disposable gloves, powderless Biohazard waste container
- Refrigerator and Freezer
- Tube racks

- Carefully read this instruction before starting the procedure.
- · For in vitro diagnostic use only.
- This assay needs to be carried out by skilled personnel.
- Clinical samples should be regarded as potentially infectious materials and
- should be prepared in a laminar flow hood.
- This assay needs to be run according to Good Laboratory Practice.
- · Do not use the kit after its expiration date.
- · Avoid repeated thawing and freezing of the reagents, this may reduce the sensitivity of the test.
- · Once the reagents have been thawed, vortex and centrifuge briefly the tubes before use.
- · Quickly prepare the reaction mix on ice or in the cooling block.
- Set up two separate working areas: 1) Isolation of the RNA/ DNA and 2) Amplification/ detection of amplification products.
- · Pipets, vials and other working materials should not circulate among working units.
- · Use always sterile pipette tips with filters.
- Wear separate coats and gloves in each area

8. Sample Collection, Storage and transportation

- Collect samples in sterile tubes;
- Specimens can be extracted immediately or frozen at -20°C to -80°C.
- · Transportation of clinical specimens must comply with local regulations for the transport of etiologic agents

9. Procedure

9.1 DNA-Extraction

DNA extraction buffer is supplied in the kit. Please thaw the buffer thoroughly and spin down briefly in the centrifuge before use

9.1.1 Stool samples

- 1) Take about 50mg samples to a 1.5ml tube; add 1.0ml normal saline then vortex vigorously. Centrifuge the tube at 13000rpm for 2 minutes, carefully remove and discard supernatant from
- the tube without disturbing the pellet.

 2) Add 100µl DNA extraction buffer, close the tube then resuspend the pellet with vortex vigorously. Spin down briefly in a table centrifuge.
- 3) Incubate the tube for 10 minutes at 100°C.
- 4) Centrifuge the tube at 13000rpm for 5 minutes. The supernatant contains the DNA extracted and can be used for PCR template

9.1.2 Water samples

- 1) Take 3 ml water to a tube, Centrifuge the tube at 13000rpm for 2 minutes, carefully remove and discard supernatant from the tube without disturbing the pellet. 2) Add $100\mu l$ DNA extraction buffer, close the tube then resuspend the pellet with vortex
- vigorously. Spin down briefly in a table centrifuge.
- 3) Incubate the tube for 10 minutes at 100°C.
 4) Centrifuge the tube at 13000rpm for 5 minutes. The supernatant contains the DNA extracted and can be used for PCR template

Attention:

- . During the incubation, make sure the tube is not open, as the vapor will
- volatilize into the air and may cause contamination if the sample is positive.
- **B.** The extraction sample should be used in 3 hours or store at -20°C for one month.
- C. Different DNA extraction kits are available. You may use your own extraction systems or the commercial kit based on the yield. For the DNA extraction, please comply with the manufacturer's instructions.


9.2 Internal Control

It is necessary to add internal control (IC) in the reaction mix. Internal Control (IC) allows the user to determine and control the possibility of PCR inhibition

Add the internal control (IC) 1µl/rxn and the result will be shown in the 560nm

9.3 PCR Protocol

The Master Mix volume for each reaction should be pipetted as follows:



※PCR system without 560nm channel may be treated with 1μl Molecular Grade Water instead of 1μl IC.

- The volumes of Reaction Mix and Enzyme Mix per reaction multiply with the number of samples, which includes the number of the controls, standards and sample prepared. Molecular Grade Water is used as the negative control. For reasons of unprecise pipetting, always add an extra virtual sample. Mix the master mix completely then spin down briefly in a centrifuge.
- 2) Pipet 18µl Master Mix with micropipets of sterile filter tips to each *Real time* PCR reaction plate/tube. Then separately add 2µl DNA sample, positive and negative controls to different reaction plate/tubes. Immediately close the plate/tubes to avoid contamination.
- 3) Spin down briefly in order to collect the Master Mix in the bottom of the reaction tubes.

refrom the following protocor in the instrument		
37°C for 2min	1cycle	
94°C for 2min	1cycle	
93°C for 5sec, 60°C for 30sec	40cycles	
(Fluorescence measured at 60°C)	,	

Selecti	Selection of fluorescence channels	
530nm	Target Nucleic Acid	
560nm	IC	

10.Threshold setting: Choose Arithmetic as back ground and none as Noise Band method, then adjust the Noise band just above the maximum level of molecular grade water, and adjust the

threshold just under the minimum of the positive control.

11.Quality control: Negative control, positive control and internal control must be performed correctly, otherwise the sample results is invalid.

v.	therwise the sample results is invalid.				
	Channel	Crossing point value			
	Control	530nm	560nm		
	Molecular Grade Water	Blank	25~33		
	Positive Control(qualitative assay)	≤35			

12. Data Analysis and Interpretation

The following results are possible:

	Crossing point value		Result Analysis	
Γ	530nm	560nm	Result Alialysis	
1#	Blank	25~33	Below the detection limit or negative	
2#	≤35		Positive	
3#	35~40	25~33	Re-test; If it is still 35~40, report as 1#	
4#	Blank	Blank	PCR Inhibition; No diagnosis can be concluded.	

