



Caspase-3 DEVD-R110 Fluorometric HTS Assay Kit

Catalog Number: 30009

Description

Caspase-3 is an active cell-death protease involved in the execution phase of apoptosis, during which cells undergo morphological changes such as DNA fragmentation, chromatin condensation, and apoptotic body formation $^{(1,2)}$. Caspase-3 DEVD-R110 Fluorometric HTS Assay Kit provides a single-step homogenous assay specifically designed for HTS-based detection. The fluorogenic substrate (Ac-DEVD)₂-R110 contains two DEVD tetrapeptides and is completely hydrolyzed by the enzyme in two successive steps. Cleavage of the first DEVD peptide results in the monopeptide Ac-DEVD-R110 intermediate, which has absorption and emission wavelengths similar to those of R110 ($\lambda_{abs}/\lambda_{abs}$ =496/520 nm), but has only about 10% the fluorescence of the latter $^{(3-4)}$. Hydrolysis of the second DEVD peptide releases the dye R110, leading to a substantial fluorescence increase.

R110-based substrate
$$\lambda_{abs}/\lambda_{em} = 232 nm/no \ emission$$

$$R110 \ h \ O \ h \ C \ C \ O_2 \ h \ C \ C \ O_2 \ h \ C \ C \ O_2 \ h \ C \ O_$$

The assay kit includes DEVD-CHO, which is a caspase-3 inhibitor and can be used as a negative control. Also, R110 is provided in the kit for generating a standard curve, which can be used for quantifying caspase-3 activity.

Kit Components

1mL (#30009-1)	10mL (#30009-2)	100mL (#30009-3)	
1mL	10mL	100mL	Cell Lysis/Assay Buffer
50uL	500uL	5mL	Enzyme Substrate (Ac-DEVD) ₂ -R110 (2mM)
5uL	20uL	100uL	Enzyme Inhibitor Ac-DEVD-CHO (5mM)
1mL	1mL	1mL	R110 (80µM)

Storage Condition

Caspase-3 DEVD-R110 Fluorometric and Colorimetric Assay Kit should be stored at –20°C or below. The components of the kit are stable at –20°C for six months. Avoid frequent freeze-thaw cycles.

Features

HTS-compatible: Single-step homogenous assay specifically designed for HTS-based detection.

Fast: Fast enzyme kinetics.

Sensitive: The enzymatic reaction forms an intensely green fluorescent rhodamine 110 (R110) product. The long wavelength of R110 excitation and emission minimize cellular autofluorescence.

Assay for Detection of Caspase-3 Activity in Cell Culture A. General Considerations

We recommend performing three control reactions:

- 1) Negative control on uninduced cells.
- 2) Control on induced cells treated with Caspase-3 inhibitor.
- 3) Positive control for Caspase-3 induction.

B. Preparation of Caspase-3 Detection Buffer

Depending on the required volume of Caspase-3 Detection Buffer, mix the Enzyme Substrate (Ac-DEVD)₂-R110 (2mM) with the Cell Lysis/Assay Buffer in a 50µL to 1mL ratio to derive Caspase-3 Detection Buffer.

C. Assay Procedure

- 1. Induce apoptosis in cells by desired methods. Remember to incubate concurrent culture without induction.
- For suspension cells, count cells and aliquot equal number of cells into each well in a 96-well plate
 or 384-well plate. It is recommended to use 500-50,000 cells per sample in the cell medium whose
 volume is equal to the volume of Caspase-3 Detection Buffer to be added. For example, cells
 should be in 100μL medium in each well if 100μL Caspase-3 Detection Buffer will be used for each
 assay.
- 3. Add Caspase-3 Detection Buffer in equal volume to cell medium directly into each well.
- 4. **[Optional]** To verify that the signal detected by the kit is due to Caspase-3 activity, incubate an induced sample with caspase-3 inhibitor before adding substrate. This can be accomplished by adding 100μL of Cell Lysis/Assay Buffer and 2μL of Enzyme Inhibitor Ac-DEVD-CHO (5mM) to the cell suspension in a well of a 96-well plate. Incubate on ice for 30 min or RT for 15 min followed by adding 5μL Enzyme Substrate (Ac-DEVD)₂-R110 (2mM).
- 5. Incubate at 37°C for 30 min to 1hr (or up to 3 hours maximum) in an incubator.
- 6. Read on a fluorometer with 470 nm excitation filter and 520 nm emission filter for optimal sensitivity.
- 7. Use R110 if necessary for generating a standard curve to calculate amount of substrate conversion.

References

- 1. Porter AG, Janicke RU. Emerging roles of caspase-3 in apoptosis. Cell Death Differ. 1999 Feb;6(2):99-104.
- 2. Zou H, Li Y, Liu X, Wang X. An APAF-1.cytochrome c multimeric complex is a functional apoptosome that activates procaspase-9. J Biol Chem. 1999 Apr 23;274(17):11549-56.
- 3. An S, Zheng Y, Bleu T. Sphingosine 1-phosphate-induced cell proliferation, survival, and related signaling events mediated by G protein-coupled receptors Edg3 and Edg5. J Biol Chem. 2000 Jan 7;275(1):288-96.
- 4. Hug H, Los M, Hirt W, Debatin KM. Rhodamine 110-linked amino acids and peptides as substrates to measure caspase activity upon apoptosis induction in intact cells. Biochemistry. 1999 Oct 19;38(42):13906-11.

Gentaur Molecular Products Voortstraat 49 1910 Kampenhout, Belgium