DATA SHEET



DNase I (RNase-Free) Deoxyribonuclease I, Ribonuclease & Protease Free

Catalog No. Quantity 15601-100 2500 unit

Description:

Endonuclease which digests dsDNA by hydrolysis of deoxyribonucleotide linkages, producing 3'-hydroxyl oligonucleotides.

Source: Bovine Pancreas

Form: Lyophilized in vials and chromatographically purified to remove RNase

and protease.

DNase Activity: ≥2,000 units/mg dry weight

Protease: None detected RNase: None detected

Unit Definition: 1 unit causes an increase in absorbance at 260nm of 0.001 per minute per ml at 25°C when acting upon highly polymerized DNA at pH 5.0. Note: 0.005 Kunitz units digests 1ug DNA in 10 minutes at 37°C, 50mM Tris,

1mM Ca+2, 1mM Mg+2 in a 50ul reaction.

Assay Method: Based on the Kunitz method where DNase will depolymerize DNA; 0.005 Kunitz units digests 1ug of lambda DNA in 10 minutes at 37°C in 50mM Tris, 1mM Ca+2 and 1mM Mg+2, pH 7.8 in a 50ul reaction. This product is for research purposes only. Not for diagnostic use.

Storage Conditions

Store at 2-8°C lyophilized for up to 2 years.

Recommended DNase I Solution Preparation and Storage
Prepare DNase I solution by adding 500 µl of sterile water to the DNase I
(RNase-Free) lyophilized powder and mix gently. Aliquotting the DNase I
stock enzyme is recommended to avoid repeated freeze/thaw cycles. The
DNase I stock enzyme can be freeze/thawed up to three times without loss of
activity.

Storage

Remove DNase I (RNase-Free) lyophilized powder and store at 4ºC. Store resuspensed DNase I solution at -20°C (DNase I is sensitive to physical denaturation. Do not vortex the resuspended DNase I solution). Store the 10x DNase I Buffer at room temperature (15-30°C) or 4°C.

DNase I Digestion Protocol for Purified RNA in Solution
This protocol is for the removal of genomic DNA from RNA samples. The reaction may be scaled up or down according to the volume of RNA in solution.

For a 100 µl reaction, prepare as follows:

RNA sample: ≤88 µl 10X buffer: 10 µl

DNase I stock enzyme: 2 µl

Sterile water: bring the volume up to 100 µl

Mix and incubate at room temperature (22-25°C) for 10 minutes.

Inactivate DNase I by adding EDTA to a final concentration of 5 mM to chelate

the divalent cations then heat to 75°C for 5 minutes.