

# Free PSA IRMA

## Catalog #: RVRK-85CT

The technology uses two high affinity monoclonal antibodies in an immunoradiometric assay (IRMA) system. The <sup>125</sup>I labelled signal-antibody binds to an epitope of the fPSA molecule spatially different from that recognized by the biotincapture- antibody. The two antibodies react simultaneously with the antigen present in standards or samples, which leads to the formation of a capture antibody - antigen - signal antibody complex, also referred to as a "sandwich". During a 2-hour incubation period with shaking immuno-complex is immobilized to the reactive surface of streptavidin coated test tubes. Reaction mixture is then discarded, test tubes washed exhaustively, and the radioactivity is measured in a gamma counter. The concentration of antigen is directly proportional to the radioactivity measured in test tubes. By constructing a calibration curve plotting binding values against a series of calibrators containing known amount of fPSA, the unknown concentration of fPSA in patient samples can determined.

<b>Format</b>	IRMA
<b>Label</b>	I125
<b>Size</b>	100 tests
<b>Sample Type</b>	Serum
<b>Sample Volume</b>	100 µL
<b>Controls</b>	1 level
<b>Range</b>	0,1-100 ng/mL
<b>Sensitivity</b>	0,047 ng/mL
<b>Incubation</b>	2 h RT with shaking
<b>Shelf Life (weeks)</b>	8
<b>Normal Values</b>	<0,7 ng/mL