

DATA SHEET

Revision No.: ZJ0001

Issue Date: Jan 18th, 2011

Clostridium perfringens Real Time PCR Kit Cat. No.: DD-0167-02

For use with ABI Prism*7000/7300/7500/7900; Smart CyclerII;Cycler iQTM4/iQTM5; Rotor GeneTM6000; Mx3000P/3005P;MJ-Option2/Chromo4;LightCycler*480real time PCR systems

User Manual For in vitro Diagnostic use only

Manufactured by Shanghai ZJ Bio-Tech Co., Ltd. www.liferiver.com.cn Tel: +86-21-34680596

trade@liferiver.com.cn Fax: +86-21-34680595 2nd floor,No.15 Building,No.188 Xinjunhuan road, PuJiang Hi-tech Park Shanghai China

A. Internet USE
Clostridium perfringens real time PCR kit is used for the detection of Clostridium perfringens in stool or water samples by using real time PCR systems.

 Principle of Real-Time PCR
The principle of Real-Time PCR
The principle of the additional time PCR.

or water samples by using real time PCR systems.

2. Principle of Real-Time PCR

The principle of the real-time detection is based on the fluorogenic 5 nuclease assay. During the PCR reaction, the DNA polymerase cleaves the probe at the 5' end and separates the reporter dye from the quencher dye only when the probe hybridizes to the target DNA. This cleavage results in the fluorescent signal generated by the cleaved reporter dye, which is monitored real-time by the PCR detection system. The PCR cycle at which an increase in the fluorescence signal is detected initially (C) is proportional to the amount of the specific PCR product. Monitoring the fluorescence intensities during Real Time allows the detection of the accumulating product without having to re-open the reaction tube after the amplification.

3. Product Description

Clostridium perfringens (formerly known as C. welchii) is a Gram-positive, rod-shaped, anaerobic, spore-forming bacterium of the genus Clostridium. C. perfringens is ubiquitous in nature and can be found as a normal component of decaying vegetation, maine sediment, the intestinal tract of humans and other vertebrates, insects, and soil. C. perfringens is a human pathogen sometimes, and other times it can be ingested and not cause any harm.

Clostridium perfringens through polymerase chain reaction in the real-time PCR system. The master contains reagents and enzymes for the specific amplification of the Clostridium perfringens DNA. Fluorescence is emitted and measured by the real time systems' optical unit during the PCR. The detection of amplified Clostridium perfringens DNA fragment is performed in fluorimeter channel FAM with the fluorescent quencher BHQ1. In addition, the kit contains a system to identify possible PCR inhubition by measuring the HEXVIC/IOE fluorescence of the internal control (C). An external positive control defined as 1×10 copies/ml is supplied which allow the determination of the gene load. For further information, please refer to section 9.3Quantitation.

4. Kit C

ontents		
Ref.	Type of Reagent	Presentation 25rxns
1	DNA Extraction Buffer	2 vials, 1.5ml
2	C. perfringens Reaction Mix	1 vial, 950μl
3	PCR Enzyme Mix	1 vial, 12μl
4	Molecular Grade Water	1 vial, 400μl
5	Internal Control	1 vial, 30μl
6	C perfringens Positive Control(1×10 ⁷ copies/ml)	1 vial 30ul

Analysis sensitivity: 1×10³ copies/ml; LOQ: 2×10³~1×10⁸ copies/ml

- storage

 All reagents should be stored at -20° C. Storage at $+4^{\circ}$ C is not recommended.

 All reagents can be used until the expiration date indicated on the kit label.

 Repeated thawing and freezing (-3x) should be avoided, as this may reduce the sensitivity of the

Real time PCR system
Real time PCR reaction tubes/plates
Pipets (0.5µl – 1000µl)
Sterile microtubes
Biohazard waste container
The nead to the property of the

- assay. Cool all reagents during the working steps.

- Reaction mix should be stored in the dark.

 Additionally Required Materials and Devices
 Biological cabinet
 Vortex mixer
 Cryo-container
 Sterile filter tips for micro pipets
 Disposable gloves, powderless
 Refrigerator and Freezer
 Desktop microcentrifuge for "eppendorf" type tubes (RCF max. 16,000 x g)
 Warnings and Precaution
 Carefully read this instruction before starting the procedure.
 For in vitro diagnostic use only.

- . Narmings and recount on ...

 Carefully read this instruction before starting the procedure.

 For in vitro diagnostic use only.

 This assay needs to be carried out by skilled personnel.

 Climical samples should be regarded as potentially infectious materials and should be prepared in a laminar flow hood.

 This assay needs to be run according to Good Laboratory Practice.

 Do not use the kit after its expiration date.

 Avoid repeated thawing and freezing of the reagents, this may reduce the sensitivity of the test.

 Once the reagents have been thawed, vortex and centrifuge briefly the tubes before use.

 Quickly repeare the reaction mix on ice or in the cooling block.

 Set up two separate working areas: 1) Isolation of the RNA/ DNA and 2) Amplification/ detection of amplification products.

 Pipets, vials and other working materials should not circulate among working units.

 Use always sterile pipete tips with filters.

 Wear separate coats and gloves in each area.

 S. Sample Collection, Storage and transportation

 Collect samples in sterile tubes;

 Specimens can be extracted immediately or frozen at -20°C to -80°C.

 Transportation of clinical specimens must comply with local regulations for the transport of enlologic agents Transportation of clinical specimens must comply with local regulations for the transport of etiologic agents
 P. Procedure
 P. Procedure
 P. I. Sampling and DNA extraction
 P. I. Sampling and increasing bacteria
 Sternlize the sample package and sampling tools before sampling
 Take 100g, 10g and 1g milk powder into three different culture bottles of 2L, 250ml and 125ml expectivable.

-) Add 9 times volume sterilized water into these three culture bottles(M/V=1:9), and incubate them under 36±1°C for 18-22h.
- 4) Take 10ml culture medium from the culture bottles into 90ml EE broth respectively, and incubate them for 18-22h.
- 9.1.2 DNA-Extraction DNA extraction buffer is contained in the kit. Please thaw the buffer thoroughly and spin down briefly in the centrifuge before use.
- in the centifing of once use.

 9.1.1 Stool samples

 1) Take about 50mg stool samples to a 1.5ml tube; add 1.0ml normal saline then vortex vigorously.

 Centrifinge the tube at 13000pm for 2 minutes, carefully remove and discard supernatant from the tube

without disturbing the pellet. 2) Add 100µl DNA extraction buffer, close the tube then resuspend the pellet with vortex vigorously.

Spin down briefly in a table centrifuge.

3) Incubate the tube for 10 minutes at 100°C.

4) Centrifuge the tube at 13000rpm for 5 minutes. The supernatant contains the DNA extracted and can be used for PCR template.

can be used for PCR template.

9.1.2 Water samples

1) Take 3 ml water to a tube. Centrifuge the tube at 13000rpm for 2 minutes, carefully remove and discard supernatant from the tube without disturbing the pellet.

2) Add 100µl DNA extraction buffer, close the tube then resuspend the pellet with vortex vigorously. Spin down briefly in a table centrifuge.

3) Incubate the tube for 10 minutes at 100°C.

4) Centrifuge the tube at 13000rpm for 5 minutes. The supernatant contains the DNA extracted and can be used for PCR template.

Attention:

Attention:

A. During the incubation, make sure the tube is not open, as the vapor will volatilize into the air and may cause contamination if the sample is positive.

B. The extraction sample should be used in 3 hours or store at -20°C for one month.

C. Different DNA extraction kits are available. You may use your own extraction systems or the commercial kit based on the yield. For the DNA extraction, please comply with the manufacturer's

instructions.

9.2 Internal Control and positive control
It is necessary to add internal control (IC) in the reaction mix. Internal Control (IC) allows the user to
determine and control the possibility of PCR inhibition.

Add the internal control (IC) | µl/rxn and the result will be got in the HEX/VIC/IOE channel.

Attention: It is necessary to dilute the internal control supplied in the kit by 10 times with molecular
grade water before detection, and close the tube immediately then vortex for 10 seconds.

Because of transportation with carbon dioxide ice, there may be white precipitate in tubes of internal
control and positive control, but it will disappear in a few minutes when it is incubated at room
temperature. Besides, the white precipitate have no effection on the detection result.


9.3 Quantitation

temperature. Besides, the white precipitate have no effection on the detection result.

9.3. Quantitation

The kit can be used for quantitative or qualitative real-time PCR. A positive control defined as 1×10² copies/mil is supplied in the kit.

For performance of quantitative real-time PCR, Standard dilutions must prepare first as follows. Molecular Grade Water is used for dilution. The step of dilution is not needed for performance of qualitative real-time PCR. Take positive control (1×10² copies/mil) as the starting high standard in the first tube. Respectively pipette 36ul of Molecular Grade Water into next three tubes. Do three dilutions as the following figures: Dilution of Standards



To generate a standard curve on the real-time system, all four dilution standards should be used and defined as standard with specification of the corresponding concentrations.

Attention:

A. Mix thoroughly before next transfer

B. The positive control (1×10⁷cc B. The positive control (1×10⁷copies/ml) contains high concentration of the target DNA.

Therefore, be careful during the dilution in order to avoid contami

The Master Mix volume for each reaction should be pipetted as follows: 1µI Internal Control 0.4µl nzyme Mix This system
is only for
Smart Cycler II PCR Ins OR

- PCR system without HEXVIC/IOE channel may be toated with jul Molecular Grade Wate unstead of jul IC

 The volumes of Reaction Mix and Enzyme Mix per reaction multiply with the number of
 samples, which includes the number of the controls. standards and sample prepared. Molecular
 Grade Water is used as the negative control. For reasons of unprecise pipetting, always add an
 extra virtual sample. Mix the master mix completely then spin down briefly in a centrifuge.

 2) Pipet 36µl (22.5µl for SmartCycer II) Master Mix with micropipets of sterile filter tips to each
 Real time PCR reaction plate/tube. Then separately add 4µl (25.µl for SmartCycer II) DNA
 sample, positive and negative controls to different reaction plate/tubes. Immediately close the
 plate/tubes to avoid contamination.

 3) Spin down briefly in order to collect the Master Mix in the bottom of the reaction tubes.

 4) Perform the following protocol in the instrument:
 37°C for 2 min, 1 cycle;94°C for 2 min, 1 cycle;93°C for 15 sec, 60°C for 60 sec, 40 cycles.
 Fluorescence is measured at 60°C;74M and HEX/TVC/JOE channels should be choren.

 5) If you use ABI Prism® system, please choose "none" as passive reference and quencher.

 10. Baseline serting; just above the maximum level of molecular grade water.

 11. Calabaration for quantitative detection: input each concentration of standard controls at the end
 of mm, and a standard curve will be automatically formed.

 12. Quality control: The Ct value of micecular grade water and positive control in FAM channel
 shows 25–33; Correlation coefficient of standard curve should be ≪-0.98, otherwise the result is
 invalid.

 3) Botto Audition and Lettomystation.

13. Data Analysis and Interpretation

The following results are possible:

1) The Ct value in channel FAM shows ≤35. The result is positive: The sample contains

1) The Ct value in channel FAM shows ≤35. The result is positive: The sample contains Clostridium perfringers DNA.

2) The Ct value in channel FAM shows 35-40, please repeat again. If the result still shows 35-40, it can be considered negative.

3) In channel FAM no signal is detected, at the same time, a HEX/VIC/JOE signal from the Internal Control appears. The sample does not contain any Clostridium perfringers DNA. It can be considered negative.

4) Neither in channel FAM nor in channel HEX/VIC/JOE signal is detected. A diagnostic statement can not be made. Inhibition of the PCR reaction.

For further questions or problems, please contact our technical support at trade@liferiver.com.cn

