



Certificate of Analysis

This is to verify that all items noted below are in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation: that all process certifications, material certifications, are on file at this or other Bel – Art facility and are available for review by the customer and or Government Representatives upon request.

Product #: 37876-0002

Description: Bel-Blotter™ 96-Well Replicating Tool

Material: Polycarbonate

Approximate Overall 5.11" x 3.42" x 1.15", Tips .65" L x 0.12" taper to 0.28" OD,

Dimensions: ID 0.074 to 0.23", Tips are spaced 0.35" center to center, volume 10ul

Other: Autoclavable at 121 C, materials meet FDA Regulation 21CFR177.1550, 21CFR177.1520

Based on the process and materials involved in the manufacturing of this product, this product is considered to be BSE/TSE protein free

Based on the information provided to us this product does not contain or come in contact with lead, latex, Mercury, Melamine, formaldehyde or other contaminants

To the best of our knowledge based on information provided to us this product does not contain or come in contact with any animal origin.