

Product name : **CIRCULATING PUMP 5L_MINUTE**

Catalog number : **DBI-01**

ISO Class **I**

DIN4050 IP33

MWC 60°C 80% H

REF	20°C H2O		l/m	50°C Silicon		Max Inm	Extern. Mm			
	l/min	mBar		mBar	Watts		D	H	W	Kg
DBI-01	5	120	2	100	20	120mm	160	215	125	3

This kind of apparatus are developed to be introduced into a recipient containing liquid, to allow to make circulate the liquid externally, to allow circulate across the serpentins and thermostatisation chambers, and stir the liquids contained on a recipient, by propulsion system as well Metallic painted housings, and all the elements in direct contact with the liquid are made of stainless steel. Equipped with flow regulation system, and provided with a double fixation system, one is through an "U" device to be fixed on the wall of the tanks, and the other one is with a screwed Ø12mm bar supplied as accessory **DXV-12**, to allow to fix on the classical laboratory bars scaffolding that is fixed on the laboratory wall. Produced two different models, with different powers and flows.

GENTAUR

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, Belgium