



CRP Antigen

Cat. # A7512H

Important Note:	Centrifuge before opening to ensure complete recovery of vial contents.
Catalog #:	A7512H
Lot #:	14F1661
Description:	CRP Antigen C Reactive Antigen (CRP) > 99% pure
Source:	Human Plasma
Format:	Purified, Liquid
Purification:	> 99% pure (SDS-PAGE)
Concentration:	2.56mg/ml (OD280nm, E0.1%=1.7)
Buffer:	20 mM Tris, pH 8.0 ± 0.1 containing 280 mM Sodium Chloride, 5 mM Calcium Chloride
Preservative:	0.1% Sodium Azide
Applications:	Specific methodologies have not been tested using this product.
Storage:	Store at 2-8°C. DO NOT FREEZE!
Inactivation:	Not Applicable
Warning:	<p>All human source materials have tested negative for HIV1, HIV2, HCV and HBV. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.</p> <p>This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1 - 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.</p>

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.