

# Calreticulin Blocking Peptide

Catalog #: 3077BP-50

## FORMULATION:

50 µg (0.5 mg/ml) in phosphate buffered saline (PBS), pH 7.2, containing 50% glycerol, 1% BSA and 0.02% thimerosal.

## STORAGE CONDITIONS:

Store at -20°C. For long-term storage, aliquot and refreeze at -70°C. Avoid repeated freeze/thaw cycles.

## APPLICATION AND USAGE:

The peptide is used for blocking the antibody activity of Calreticulin (Cat.# 3077-100). It usually blocks the antibody activity completely in Western blot analysis by incubating the peptide with equal volume of antibody for 30-60 minutes at 37°C.

***FOR RESEARCH USE ONLY! Not to be used in humans.***