



CappTronic Capp electronic pipette, 8-channel

The CappTronic™ is certainly one of the easiest electronic pipettes to set and use. Having only a one-side programming manual, familiarizing ones self with CappTronic's™ functions is a mater of minutes. Apart from standard pipetting, you can choose between the most widely used pipetting modes including reverse pipetting, multi dispensing, mixing and sequential dispensing. Up to nine individual protocols can be stored thanks to the extensive memory capacity. Supplied with two long-life and easily exchangeable lithium batteries, the CappTronic™ allows for uninterrupted pipetting without the need to stop work in order to recharge.

The O-rings on the tip bases provide an excellent fit with a great variety of universal fitting tips.

Cat No	Description	Inaccuracy %	Imprecision %	Color
E20-8BA	Capp electronic pipette, 8-channel, 2-20 ul	1.00	0.30	red 🌑
E200-8BA	Capp electronic pipette, 8-channel, 10-200 ul	0.60	0.15	yellow 🔵
E1200-8BA	Capp electronic pipette, 8-channel, 100-1200 ul	0.50	0.15	blue





Gentaur Molecular Products Voortstraat 49 1910 Kampenhout, Belgium