

Carbogen: Rapid Plasma Reagin (RPR) card test / carbon antigen for syphilis testing

Intended Use:

Carbogen, a Rapid Plasma Reagin (RPR) card test / carbon antigen for the detection of Anti-Lipoidal antibodies in serum and plasma.

Carbogen kit components:

Reagent	Microfine particulate carbon suspension coated with lipid complexes	
Positive Control	Reactive with the reagent	
Negative Control	Non Reactive with the reagent	
Slide	Disposable slide with eight reaction circles	
Other Accessories	Other Accessories Dropper, pipette, mixing sticks, rubber teats and Package Insert	

Performance Characteristics:

Sensitivity: 100 % Specificity: Lipid complex reactive with Anti-Lipoid	dal antibodies
---	----------------

FEATURES	BENEFITS	
Specially designed anti – complement agent removes complement interference	No false reactions due to complements	
Homogenous and micronised antigen carrier	Clear - cut negatives	
Coated disposable cards provides optimum reaction environment	Unmatched readability	
Specially purified cardiolipin antigen approved by CDC, USA	High sensitivity and specificity	
Reagents standardized against WHO IRP	Reliable performance of kit	
Kits exported to 45 countries worldwide	World – class quality	
Special pack design	Saves refrigerator space	

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	18 Months
Opened kit	2-8°C	18 Months

Available Pack Sizes		
5 ml	50 Tests	
100 Tests	250 Tests	