

DATA SHEET

Product Description

Product: Chicken plasmin

Catalog Number: CKPLM

Description: Prepared from plasminogen by activation with immobilized human uPA. 100 percent functionally

active plasmin is purified from the activation reaction by immobilized SBTI. Plasmin will undergo rapid auto proteolysis in the absence of benzamidine and should be used quickly once thawed.

Aliquot to avoid freeze-thaw cycles.

Recommended Storage: -70 C **Shipping Conditions:** Dry ice

Minimum Shelf Life: 3 years from delivery

Technical Information

Purity: >95% by SDS-PAGE analysis

Source: Chicken plasma

Molecular Weight: 85,000 **Extinction Coefficient:** 1.69

Sample Characteristics

Lot Number: CKPLM-SAMPLE

 Aliquot:
 1 x 1.0 mg

 Concentration:
 5.0 mg/ml

 Volume:
 0.2 ml

Buffer: 0.1M HEPES; 0.1M NaCl; pH 7.4

Form: Frozen liquid

Optimal experimental conditions and working dilutions must be determined by the end user. Sample concentration and volume are for example purposes only and do not necessarily represent the product that will be received. Please contact us for a current lot-specific product datasheet.

For in vitro laboratory use only

Gentaur Molecular Products Voortstraat 49 1910 Kampenhout, Belgium