

DATA SHEET

Product Description

Product: Chicken plasminogen

Catalog Number: CKPG

Description: Prepared from fresh chicken plasma by immobilized lysine chromatography.

Recommended Storage: -70 C **Shipping Conditions:** Dry ice

Minimum Shelf Life: 3 years from delivery

Technical Information

Purity: >95% by SDS-PAGE analysis

Source: Chicken plasma

Molecular Weight:92,000Extinction Coefficient:1.69

Sample Characteristics

Lot Number: CKPG-SAMPLE

Aliquot: $1 \times 1.0 \text{ mg}$ Concentration: 5.0 mg/mlVolume: 0.2 ml

Buffer: 0.1M HEPES; 0.1M NaCl; pH 7.4

Form: Frozen liquid

Optimal experimental conditions and working dilutions must be determined by the end user. Sample concentration and volume are for example purposes only and do not necessarily represent the product that will be received. Please contact us for a current lot-specific product datasheet.

For in vitro laboratory use only

Gentaur Molecular Products Voortstraat 49 1910 Kampenhout, Belgium