



Neptune Barrier Tips

NEPTUNE™



Protect your work from contamination with one of the industry's most efficient aerosol filter assemblies. Our proven range of filter tips reliably prevent aerosol transfer to give you complete peace of mind when performing your most sensitive assays.

The exclusive S3 polymer was designed to increase pipetting accuracy by virtually eliminating tip retention and sample hold-up. S3 pipette tips deliver volumes within 0.1% of the indicated volume, versus 0.7% for standard polymer tips. This provides a more accurate application.



Cat. Nr.	Product	Quantity	Price
BT10	NEPTUNE Barrier Tips 10µl	10x96	€ 59



BT10XL	NEPTUNE Barrier Tips 10µl XL	10x96	€ 59
--------	------------------------------	-------	------



BT20	NEPTUNE Barrier Tips 20µl	10x96	€ 59
------	---------------------------	-------	------



BT100	NEPTUNE Barrier Tips 100µl	10x96	€ 59
-------	----------------------------	-------	------



BT200	NEPTUNE Barrier Tips 200µl	10x96	€ 59
-------	----------------------------	-------	------



BT1250	NEPTUNE Barrier Tips 1250µl	8x60	€ 34
--------	-----------------------------	------	------

OrTube centrifuge tubes



1 box - ~~€35~~ €25
10 boxes - ~~€300~~ €250
20 boxes - ~~€540~~ €450
30 boxes - ~~€760~~ €600



Cat. Nr.	Product
04-5540400	15ml TC Tubes, Round, PP, 560 pcs/box
04-5540150	50ml TC Tubes, Conical, 440 pcs/box

Price
€25/box
€25/box

OrPlate culture plates



- Made of USP VI class virgin Polystyrene, high clarity and toughness, high light transferability with minimal light distortion
- Guaranteed even surface refinement in the growth zone
- Alpha-numeric identification method to identify each well
- Specially designed lid separates each well, thus preventing liquid or condensate cross-contamination
- Deep spacing between the wells to prevent cross-contamination
- Cut corner makes the lid a "one-way" lid
- Special design allows easy and stable stacking
- Optimum optical clarity for microscope observation
- Packaged in easy to peel paper back and plastic packs
- Test Plates are certified as having been sterilised by gamma radiation
- Test plates are certified as being DNA and RNA free
- Test plates are certified as having an endotoxin level of <0.5 EU/ml

Cat. nr.	Product	Pcs/bag	Bags/box	Price	Price (10 boxes)
5530300	TC Test Plate, 24 wells (85 Units/box)	1	85	€35	€235
5530100	TC Orplate Test Plate, 96 Flat, 1 Unit/bag, 115 Units/box	1	115	€35	€235



www.gentaur-worldwide.com
info@gentaur.com
Tel.: 0032 16 58 90 45

