

## **DATA SHEET**

Colchicine

H<sub>3</sub>C — C — HN — OCH<sub>3</sub>

Synonym: N-[(7S)-5,6,7,9-Tetrahydro-1,2,3,10-tetramethoxy-9-

oxobenzo[a]heptalen-7-yl]acetamide

CAS: 64-86-8Formula:  $C_{22}H_{25}NO_6$ Molecular Wt: 399.44

**Properties** 

Form: Powder

Appearance: Off-white to Yellow Powder

Application: Mutagen Solubility: Water

Typical Working Concentration:

Storage Temp: Room Temperature

Storage Temp of

Stock Solution:

Other Notes: Plant Tissue Culture Tested

## **Application Notes**

This product is commonly used to double chromosomes

## References

Merck 13, 2496

Revised 3/2007

Gentaur Molecular Products Voortstraat 49 1910 Kampenhout, Belgium