

Combi-Syringe Tip, 50.0ml
Catalog #: 33-428S

- Fits Eppendorf Repeater
- 100 Sterile Syringe Tips/Unit

Replacement syringe tips designed for Eppendorf Repeater™ and Labpette R™ repeating pipettes. Labnet Biofree Combi-Syringe Tips provide accurate, reliable pipetting when used with Eppendorf Repeater or Labpette R repeating pipettes. Manufactured from virgin polypropylene with a polyethylene inner plunger, each reservoir combi tip is graduated to indicate the total capacity and the repetitive volumes it can dispense. Volumes from 1ul to 5ml can be dispensed, depending upon the reservoir tip selected.

These sterile tips are supplied in sturdy, stackable dispenser boxes of 100 tip