

### DNM-9606 Microplate reader

#### Usage:

The instrument is mainly used to immunity detection and microanalysis in chemistry and biochemistry in hospital and relevant clinical laboratory.



#### Feature:

1. With the function of automatic self-checking, filters self-choosing, board quivering and communication.
2. Storing 200 detection item, 200 calibrating parameter and its curve; it also restores 105 test data for the same board, Among of which;
3. Support the function of single or double wavelength for users' choice.
4. With the function of absorbance, determine the nature and quantitative analysis;
5. 640\*480 color touch screen, with printer outside or inside.