

Diamine oxidase

Cat. # 900-53-0

Substance name	Diamine oxidase from porcine kidney		
	≥0.05 unit/mg solid		
Art No.	140194		
CAS No.	9001-53-0		
Storage temp.	-20°C		
Test items	Specification		
Enzymatic activity	> or = 0.05 units/mg solid		
Unit definition	One unit will oxidize 1.0 μmole of putrescine per hour at pH 7.2 at 37 °C.		
Monoamine oxidase	< or = 1%		
	April,17,2012/CH28450		
Above information do not relieve the purchaser from controlling the product on delivery.			
This certificate is an electronic copy of the certificate available in the laboratory and does not require signature			