



## DATA SHEET

---

### Product Description

<b>Product:</b>	Dog plasminogen
<b>Catalog Number:</b>	DPLG
<b>Description:</b>	Prepared from fresh US origin dog (canine) plasma by immobilized lysine chromatography.
<b>Recommended Storage:</b>	-70 C
<b>Shipping Conditions:</b>	Dry ice
<b>Minimum Shelf Life:</b>	3 years from delivery

### Technical Information

<b>Purity:</b>	>95% by SDS-PAGE analysis
<b>Source:</b>	Dog plasma
<b>Molecular Weight:</b>	92,000
<b>Extinction Coefficient:</b>	1.69

### Sample Characteristics

<b>Lot Number:</b>	DPLG-SAMPLE
<b>Aliquot:</b>	1 x 1.0 mg
<b>Concentration:</b>	10.0 mg/ml
<b>Volume:</b>	0.1 ml
<b>Buffer:</b>	0.1M HEPES; 0.1M NaCl; pH 7.4
<b>Form:</b>	Frozen liquid

Optimal experimental conditions and working dilutions must be determined by the end user. Sample concentration and volume are for example purposes only and do not necessarily represent the product that will be received. Please contact us for a current lot-specific product datasheet.

For *in vitro* laboratory use only

Gentaur Molecular Products  
Voortstraat 49  
1910 Kampenhout, Belgium