



DATA SHEET

Product Description

Product:	Donkey IgG, Protein G Purified
Catalog Number:	DK-GF
Description:	Purified from normal serum by immobilized Protein G. >95 percent pure by SDS-PAGE and preservative free.
Recommended Storage:	-70 C
Shipping Conditions:	Dry ice
Minimum Shelf Life:	5 years from delivery

Technical Information

Purity:	IgG fraction
Source:	Donkey serum
Molecular Weight:	160,000
Extinction Coefficient:	1.36
Target:	This antibody has been selected as an isotype control as it has no known specificity.
Host:	Donkey
Clonality:	Polyclonal
Isotype:	IgG
Applications:	ELISA, WB

Sample Characteristics

Lot Number:	DK-GF-SAMPLE
Aliquot:	1 x 50 mg
Concentration:	10 mg/ml
Volume:	5.0 ml
Buffer:	0.02M Sodium Phosphate; 0.15M NaCl; pH 7.4
Form:	Frozen liquid

Optimal experimental conditions and working dilutions must be determined by the end user. Sample concentration and volume are for example purposes only and do not necessarily represent the product that will be received. Please contact us for a current lot-specific product datasheet.

For *in vitro* laboratory use only

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, Belgium