



DATA SHEET

Product Description

Product:	Ecotin (E. coli serine protease inhibitor)
Catalog Number:	ECOTIN
Description:	Serine protease inhibitor purified from E. coli. At least 98percent pure as determined by SDS-PAGE. References: 1. Seymour, J.L. et al. (1994) Biochemistry 33:3949-3958.
Recommended Storage:	-70 C
Shipping Conditions:	Dry ice
Minimum Shelf Life:	3 years from delivery

Technical Information

Purity:	98% by SDS-PAGE analysis
Source:	E. coli
Molecular Weight:	16,096 (monomer)
Extinction Coefficient:	1.37

Sample Characteristics

Lot Number:	ECOTIN-SAMPLE
Aliquot:	1 x 1.0 mg
Concentration:	2.5 mg/ml
Volume:	0.4 ml
Buffer:	1mM HCl
Form:	Frozen liquid

Optimal experimental conditions and working dilutions must be determined by the end user. Sample concentration and volume are for example purposes only and do not necessarily represent the product that will be received. Please contact us for a current lot-specific product datasheet.

For *in vitro* laboratory use only

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, Belgium