

# Apolipoprotein A-1, 25-267 aa, Human

Catalog Number: AS03632

**Description** Apolipoprotein A-1, 25-267 aa, Human, His-tagged, Recombinant, *E.coli*

**Form** Liquid. In 20 mM Tris-HCl buffer (pH8.0) containing, 10% glycerol

**Size** 0.1 mg

**Preparation** Recombinant, *E. coli*

**Concentration** 1.1 mg/ml (0.091 ml)


**Sequence** MGSSHHHHHH SSGLVPRGSH MDEPPQSPWD RVKDLATVYV DVLKDSGRDY  
VSQFEGSALG KQLNLKLLDN WDSVTSTFSK LREQLGPVTQ EFWDNLEKET EGLRQEMSKD  
LEEVKAKVQP YLDDFQKKWQ EEMELYRQKV EPLRAELQEG ARQKLHELQE KLSPLGEEMR  
DRARAHVDAL RTHLAPYSDE LRQLAARLE ALKENG GARL AEYHAKATEH LSTLSEKAKP ALEDLRQGLL  
PVLESFKVSF LSALEEYTKK LNTQ

**Molecular Weight** 30.3 kDa (264aa), confirmed by MALDI-TOF

**Purity** > 95% by SDS-PAGE

**Storage** Store at 4 °C short term (1-2 weeks) Store at <-20°C long term

**Gel Image**



FOR RESEARCH USE ONLY

---