

MANNOSE RECEPTOR. CD206, Antibody

GenWay ID:	GWB-57D88E
Legacy ID:	20-321-175139
Size:	0.1 mg
Source:	Mouse IgG1
Reactivity:	Human

Details for MANNOSE RECEPTOR. CD206, Antibody

Clone: 15-2 **Applications:** The antibody can be used for immunohistology on frozen sections. immuno assays and immuno precipitation. Furthermore the antibody is useful for flow cytometry. Western blotting and ligand-receptor intervention studies. The Mannose Receptor (MR. also

known as CD206). a member of the vertebrate C-type lectin family. is a pattern recognition receptor that is involved in both innate and adaptive immunity. The 180 kDa transmembrane protein consists of 5 domains: an amino-terminal cysteine-rich region. a fibronectin type II repeat. a series of eight tandem lectin-like carbohydrate recognition domains (responsible for the recognition of mannose and fucose). a transmembrane domain. and an intracellular carboxy-terminal tail. The structure is shared by the family of multi lectin mannose receptors: the phospholipase A2-receptor. DEC 205 and the novel C-type lectin receptor (mannose receptor X). The MR recognises a wide range of gram positive and gram negative bacteria. yeasts. parasites and mycobacteria. The MR has also been shown to bind and internalize tissue-type plasminogen activator. MR's are present on monocytes and dendritic cells (DC) and are presumed to play a role in innate and adaptive immunity. the latter via processing by DC. The expression of MR as observed in immunohistology is present on tissue macrophages. dendritic cells. a subpopulation of endothelial cells. Kupffer cells and sperm cells. The expression of MR on monocytes increases during culture and can be enhanced by cytokines as IFN-gamma. Labeling of MR expressing monocytes/macrophages increases with prolonged incubation time probably due to internalization of the MR-antibody-complex. The antibody 15-2 prevents binding of glycoproteins including t-PA to MR. Detection of the MR with anti-MR monoclonal antibody 15-2 can substitute staining for mannose containing probes as labeled mannosylated BSA. a technique which is more cumbersome and less specific.

N/A **[1]** Fanibunda,S.E., Bandivdekar,A.H., Velhal,S.M. and Raghavan,V.P., et al.
Mannose receptor as a CD4 independent HIV binding protein on human spermatozoa

[2] Fanibunda,S.E., Bandivdekar,A.H., Velhal,S.M. and Raghavan,V.P., et al.
Direct Submission

[3] Fanibunda S.E., Bandivdekar A.H., Velhal S.M., Raghavan V.P.
Mannose receptor as a CD4 independent HIV binding protein on human spermatozoa.

Additional Info for MANNOSE RECEPTOR. CD206, Antibody

Related Product Names	MONOCLONAL ANTIBODY TO HUMAN MANNOSE RECEPTOR. CD206 MANNOSE RECEPTOR. CD206
NCBI Acc Number	ABG47462.1
GI Number	109895388
Format	1 ml (100 ug/ml) 0.2 um filtered antibody solution in PBS. containing 0.1% bovine serum albumin.
Stability	Stable for 1 year
Clonality	Monoclonal
Storage	Product should be stored at 4 degree C
Molecular Weight	165978
Sequence Length	1456

Swiss Prot Number

Q0Z8D6