

# Mouse Thrombocyte (Platelet), Adsorbed, Rabbit anti- Catalog #: AIAD31440



<u>Product #</u>	AIAD31440
	Adsorbed Rabbit anti-Mouse Thrombocyte
<u>Test Type:</u>	For Cytotoxic Antibodies: Modified Colorimetric Microtiter Assay (Green et.al., J. Imm. Methods, vol. 70: 257, (1984))
	For Agglutinating Antibodies: Antisera dilutions in RPMI-1640 incubated with target cells at 4-8°C for 1 hr. Agglutination determined by microscopic observations.
<u>Results:</u>	Antisera of the above lot strongly agglutinates <u>mouse thrombocytes but not thymocytes or spleen cells</u> at dilutions: <u>1:1000.</u>
<u>Sterility tests:</u>	0.2 um filtered
<u>Preservatives:</u>	none
<u>Storage Conditions:</u>	Divide into useful aliquots and store at -20°C. Avoid repeated freeze/thaw steps.

