



Datasheet: **ABS3307**

Description:	RABBIT ANTI HUMAN C4d
Specificity:	C4d
Other names:	COMPLEMENT COMPONENT 4d
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.25 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	■			1/20 - 1/80
Immunohistology - Paraffin(1)	■			1/20 - 1/80
Immunofluorescence	■			1/20 - 1/80

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

(1) Heat induced epitope retrieval is recommended.

Target Species	Human
Product Form	Purified Ig - liquid
Preparation	Purified Ig prepared by affinity chromatography on Protein G
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	0.2 mg/ml
Immunogen	Synthetic peptide corresponding to amino acids 1241-1256 of C4 conjugated to Keyhole limpet Haemocyanin.
External Database Links	UniProt: P0C0L4 CO4A_HUMAN Related reagents P0C0L5 CO4B_HUMAN Related reagents Entrez Gene: 720 C4A Related reagents 721 C4B Related reagents
Synonyms	CO4, CPAMD2, CPAMD3

References	1. Regele, H. <i>et al.</i> (2001) Endothelial C4d deposition is associated with inferior kidney allograft outcome independently of cellular rejection. Nephrol. Dial. Transplant 16: 2058-2066.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet Documentation #10040 available at: http://www.abdserotec.com/uploads/MSDS/10040.pdf

GENTAUR