



## Anti Mouse Pou1F1 (dw; Hmp1; Pit1; GHF-1; Pit-1; dwarf; Pit1-rs1; Pou1f1)

<b>Product Name</b>	Rabbit Polyclonal Anti Mouse : Pou1 F1
<b>Price</b>	150 € / 500ul
<b>Antigen</b>	Séquence Peptidique d'environ 60 acide amine de transporteur Na/Pi
<b>Primary source</b>	MGI:97588
<b>Gène type</b>	Protéin Coding
<b>See related</b>	Ensembl:ENSMUSG00000004842
<b>Usage</b>	Western Blot ( 1/500 ) , ELISA, immunoprecipitation.
<b>Raised in</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Form</b>	Crude serum
<b>Storage instructions</b>	Aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
<b>Expiration</b>	20/04/2016
<b>Reference</b>	Christine Leroy Centre de Recherche « Croissance et Signalisation – Paris »

**Official Symbol** Pou1f1

**Official Full Name** POU domain, class 1, transcription factor 1

**Organism** Mus musculus

**Lineage** Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Sciurognathi; Muroidea; Muridae; Murinae; Mus

**Also known as** dw; Hmp1; Pit1; GHF-1; Pit-1; dwarf; Pit1-rs1; Pou1f1

