



DATA SHEET

Product Description

Product:	Rabbit anti human factor X antiserum
Catalog Number:	ASHFX
Description:	Polyclonal antiserum (host rabbit) to human factor X.
Recommended Storage:	-70 C
Shipping Conditions:	Dry ice
Minimum Shelf Life:	5 years from delivery

Technical Information

Purity:	Serum
Source:	Rabbit serum
Antigen:	Full length native protein (purified) (Human).
Target:	Recognises Factor X
Target Molecular Weight:	58,900
Host:	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Applications:	ELISA, WB
Cross-reactivity:	Hu

Sample Characteristics

Lot Number:	ASHFX-SAMPLE
Aliquot:	1 x 1.0 ml
Concentration:	Whole Antiserum
Volume:	1.0 ml
Form:	Frozen liquid

For In Vitro laboratory use only

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, Belgium