



VEGF-C Polyclonal Antibody

VEGF-C Antibody

Catalog#: 5633-100

Size: 100 µg

Antibody Target: VEGF-C

Alternative Name: VEGF C, Vascular Endothelial Growth Factor C , Vascular endothelial growth factor-related protein C; VRP

Host: Goat

Antibody Type: Polyclonal

Isotype: Goat IgG

Immunogen: Recombinant human VEGF-C

Accession#: P49767

Gene ID: 7424

Appearance: Colorless liquid

Concentration: 0.5 mg/ml

Formulation: 100 µg (0.5 mg/ml) antigen affinity purified goat polyclonal antibody in phosphate-buffered saline (PBS) containing 50% glycerol, 1% BSA, and 0.02% thimerosal.

Purification: Affinity purified

Species Reactivity: human

Application: Western blot, Immunohistochemistry

Application & Usage: Western blotting (0.5-1 µg/ml) and Immunohistochemistry with both frozen and paraffin sections (15-25 µg/ml). However, the optimal concentrations should be determined individually. The antibody recognizes human VEGF-C. Reactivity to other species has not been tested.

Storage Temp.: -20 °C

Shipping: gel pack

Background Descriptions:

Handling: The antibody solution should be gently mixed before use.

Usage: For Research Use Only! Not to be used in humans.

| | | |
|--------------------------|----------------------------|-------------------------------------|
| Related Products: | Secondary Antibodies | Phospho-Specific Antibodies |
| | Blocking Peptides | Protein A/G/L & Affinity Conjugates |
| | Loading Control Antibodies | Tag Antibodies |
| | Marker Antibodies | Western Blot Control Lysates |
| | Mounting Media | Western Blot Substrate Kit |
| | | |

GENTIA