

Rabbit Polyclonal to SYTL4.

Category	Primary antibodies
Product name	SYTL4 antibody
Catalog Number	orb128809
Host	Rabbit
Clonality	Polyclonal
Tested applications	WB, P-ELISA
Concentration	1mg/ml
Reacts with	Human, Mouse, Rat
Conjugation	Unconjugated
Purity	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Form	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Aliquot and store at -20 °C. Avoid repeated freeze/thaw cycles.
Note	For research use only.

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, BELGIUM

Tel 0032 16 58 90 45 | Fax 0032 16 50 90 45
www.gentaur-worldwide.com
info@gentaur.com