

Rabbit Polyclonal to TJP2.

Category	Primary antibodies
Product name	TJP2 antibody
Catalog Number	orb140339
Host	Rabbit
Clonality	Polyclonal
Tested applications	WB
Dilution range	WB: 1:1000
Reacts with	Human
NCBI	NP_004808.2
Conjugation	Unconjugated
Immunogen	This TJP2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 980-1009 amino acids from the C-terminal region of human TJP2.
Purity	Peptide Affinity Purified.
MW	133972 Da
Form	Supplied in PBS with 0.09% (W/V) sodium azide.
Application notes	Gene Symbol: TJP2. Gene-IDs: 9414. Accessions: Q9UDY2, NP_004808.2, NP_001164101.1.
Storage	TJP2 Antibody (C-term) can be refrigerated at 2-8°C for up to 6 months. For long term storage, place the at -20°C.
Note	For research use only.

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, BELGIUM

Tel 0032 16 58 90 45 | Fax 0032 16 50 90 45
www.gentaur-worldwide.com
info@gentaur.com