

Peptide affinity purified rabbit polyclonal antibody.

Category	Primary antibodies
Product name	CLTB antibody
Catalog Number	orb30957
Host	Rabbit
Clonality	Polyclonal
Tested applications	WB
Dilution range	WB: 1:1000
Reacts with	Human
NCBI	NP_009028.1
Conjugation	Unconjugated
Immunogen	Immunogen is CLTB.
Purity	Peptide affinity purified.
Isotype	IgG
Form	Supplied in PBS with 0.09% (W/V) sodium azide.
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only.

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, BELGIUM

Tel 0032 16 58 90 45 | Fax 0032 16 50 90 45
www.gentaur-worldwide.com
info@gentaur.com